

### Application Sheet 89 SPECTROFLAME-ICP

## Wear metals in used engine oils



### Introduction

Periodic engine oil analysis is used to determine which engine components are wearing and if the level of wear is nearing a critical point: drain intervals are extended and engine life prolonged.

### Application

The detection limits for Al, Cd, Cr, Cu, Fe, Mg, Mn, Mo, Na, Ni, P, Pb, Si, Sn, Ti and V are measured with ICP-OES after dilution with kerosene.

### Instrumentation

All measurements are obtained using a SPECTROFLAME equipped with an ICP 27 MHz free running source unit. The selected operating conditions are listed in Table 1.

**Table 1 - Operating conditions**

Power	1.2 kW
Gas flows	
• cooling	15 l/min
• auxiliary	1.5 l/min
Torch	
• fixed	
• Injector tube 2 mm	
Nebulizer	
• type	cross-flow
• pressure	3 bar
• carrier	0.8 l/min
• uptake rate	1.5 ml/min

### Sample preparation

Dilution 1:9 reduces viscosity effects between the different oils samples and allows an efficient introduction.

### Results

Calibration lines are established for each element using Conostan S21. Regression analysis is carried out to evaluate the detection limits.

Results are shown in Table 2.

**Table 2 - Line selection and detection limit**

Element	Wavelength (nm)	Detection limit $\mu\text{g/g}$ in oil
Al	308.21	0.1
Cd	226.50	0.04
Cr	267.72	0.09
Cu	324.75	0.1
Fe	259.94	0.07
Mg	279.55	0.003
Mn	257.61	0.01
Mo	281.61	0.2
Na	588.99	2.2
Ni	231.60	0.2
Pb	220.35	0.4
Si	251.61	0.3
Sn	189.99	0.4
Ti	337.27	0.06
V	311.07	0.09

